



THE EASY,
DURABLE WAY TO STAY DRY®

Seal-Tight Sport Instructions for Use

Did you know you can keep your cast or bandage dry while participating in your favorite water activity? You can with Seal-Tight Sport....The easy durable way to stay dry.

Here's how easy it is!

Instructions for Use

1. After removing your cast protector from the box, unfold completely.
2. Get ready for your new water activity by stretching the opening of the cast protector over your cast while pulling up. You will know it's at the right position on your arm or leg when the seal opening sits comfortably above the top cast edge.

Important note: The seal should be tight, after all, it's keeping water off of your cast. Never trim with scissors or puncture. It's also ok if the bag appears oversized, as this just allows room for large casts and easy movement.

3. Being careful not to submerge your leg past the seal, insert protector into the water while simultaneously holding a portion of the seal away from your skin with one finger. This will allow any excess air to escape easily from the bag.
4. When the seal is approximately one inch from the water's edge, remove your finger so the seal fits snugly against your skin and you're ready to go.



Removal Instructions

After swimming, removal is easy!

Just dry the protector thoroughly before stretching the seal over cast and remove.

Used properly, your Seal-Tight Sport will last the life of your cast or dressing, help prevent infection, and remove the need for costly cast replacement.

Always check with your doctor to be sure swimming is appropriate in your situation and enjoy your Seal-Tight Sport.



Scan for
product demo



www.brownmed.com

Executive:

280 Summer Street
Suite 400
Boston, MA 02210
P 857-317-3354
F 857-317-3355

Sales:

4435 Main Street
Suite 820
Kansas City, MO 64111
TF 877-853-5518
P 816-581-7001
F 816-581-7005

Manufacturing:

1300 Lundberg Drive West
Spirit Lake, IA 51360
TF 800-843-4395
P 712-336-4395
F 712-336-2874

©2010 Brownmed - SEAL-TIGHT, the SEAL-TIGHT logo, "The easy durable way to stay dry" and "Stay dry while getting wet" are trademarks or registered trademarks of Brownmed.