



Vibrating Step Sensor™



USE FOR

Gait abnormality associated with peripheral neuropathy and other conditions

WHY?

May increase walking speed, improve steadiness and increase confidence

Pressure switch applies easily to insole of any athletic or dress shoe


Provides a gentle reminder in the form of vibration as user takes a step

Hook and loop closure

Recharges easily



Scan for a product demo or visit brownmed.com/videos.



www.brownmed.com



Executive:

280 Summer Street
Suite 400
Boston, MA 02210
P 857-317-3354
F 857-317-3355

Sales:

4435 Main Street
Suite 820
Kansas City, MO 64111
TF 877-853-5518
P 816-581-7001
F 816-581-7005

Manufacturing:

1300 Lundberg Drive West
Spirit Lake, IA 51360
TF 800-843-4395
P 712-336-4395
F 712-336-2874

Item #	Description
07105	One size - Universal
RM07502	Charger